



GerMedUSA Inc.
THE COMPLETE SHOP FOR SURGICAL INSTRUMENTS




Why should you buy **from us?**

because, we not only sell, but also

manufacture Surgical Instruments

hoofcare instruments



2417 Jericho Turnpike, #333, Garden City Park, NY 11040.
Tel: 800-330-1322 - 516-358-2180 Fax: 800-260-7177 - 516-358-2182
email: sales@germdusa.com, www.GerMedUSA.com



Hoof Nippers

- V90-1001 Nipper 15" \$149.00
- V90-1002 Nipper 14" \$156.00
- V90-1004 Nipper 12" \$139.00



Carbide Inserted Jaws for Durability



Hoof Pull Offs

- D90-1200 Pull Off 14" \$129.00
- D90-1201 Pull Off 12" \$129.00



Carbide Inserted Jaws for Durability



Half-Round Nipper

- V90-1200 Half Round Nipper 12" ...\$179.00



Nail Puller

- V90-1316 Hoof Crease nail puller . \$69.00



Hoof Testers

- V90-1100 Large Hoof tester 18" \$165.00
- V90-1101 Small Pony Hoof Tester 15".... \$165.00



Hoof Clinchers

- V90-1300 Nail Curved Clincher \$130.00
- V90-1301 Race Track Clincher \$130.00



All our instruments are made of Surgical grade stainless steel. Hoof Nippers, Pull Offs have Carbide inserted jaws which makes the metal 17 times stronger than the Surgical grade Stainless Steel.



Hoof Searcher / Hoof Groover

V90-1310 Hoof Groover \$29.00
V90-1312 Hoof Searcher \$29.00



Hoof Picker

V90-1309 Hoof Picker small \$15.50



Hoof Buffer

V90-1315A Hoof Buffer \$42.00



Hoof Parer

G1107-02 Hoof Parer 14" \$119



Hoof Rasp

G1106-01 14"x1 3/4"\$43.36



Hoof Knife Kit

V90-1317 Hoof Knife Kit \$49.00



Foot-Rot Hoof Shear

V90-1315 Foot-Rot Hoof Shear s/s ... \$21.60



Hoof Knife / Blade

V90-1200B Folding Hoof Knife \$19.97
V90-1308 Swiss Hoof Knife\$25.50



Hoof Knives




Right handed



Left handed

V90-1200
Hoof Knife Regular
 Regular blade. This durable blade is held by two rivets in the hardwood handle and it extends 2" into handle. 1/2" blade at heel, 3" long. 8.5" long overall. Wooden grip.



Right handed




Left handed

V90-1200-1
Hoof Knife
 Double edged. Polished wooden grip.



Right handed




Left handed

V90-1200-2
Hoof Knife
 Polished wooden grip.



Right handed




Left handed

V90-1200-3
Hoof Knife
 Narrow blade, Polished wooden grip.



Right handed



Left handed

V90-1200-4
Hoof Knife
 Double edged, Polished wooden grip.



Right handed



Left handed

V90-1200-5
Hoof Knife
 Curved blade, Polished wooden grip.



Right handed



Left handed

V90-1200-10
Hoof Knife
 Extra long shaft makes it easier and more accurate to use. Polished Wooden grip.



Right handed




Left handed

V90-1200-11
Hoof Knife
 Narrow blade, Fine Wooden grip.



Right handed



Left handed

V90-1200-12
Hoof Knife
 Wide blade, Fine Wooden grip.



V90-1200-13
Hoof Knife
 Double edged. Length 7.5" s/s. polished Wooden grip.



V90-1200-14
Hoof Knife Small Size
 Double edged. Length 6.5" s/s. polished Wooden grip.



V90-1311
Oval Hoof Knife
 Available in wooden handle, can be used left or right handed.



V90-1311A
Oval Hoof Knife
 Available in wooden handle, can be used left or right handed. Extra long shaft makes it easier and more accurate to use.